



1. Touchless probe cover loading and ejection supports infection control.
2. Probe design resists breakage and increases durability.
3. Multiple measurement site and mode options (oral, axillary, rectal, direct, and cold) offer flexibility to obtain the most accurate and clinically relevant measurements.
4. Icon based user interface on backlit LCD display provides ease of use.

Accurate and Fast

The FILAC™ 3000 Electronic Thermometer product family offers best in class accuracy and speed. Our multiple thermometer formats address the need for accurate, fast reading in under 6 seconds. Our patented system design, controlled calibration methods and stringent manufacturing process contribute to our outstanding accuracy and reliability.

Easy to Use

The FILAC™ 3000 Electronic Thermometer's compact, robust design will stand upright on a flat surface and rubber side grips make it easy to hold. The intuitive icon based interface system with backlit LCD display is easy to read and our unique progress bar indicates status temperature measurement. AD model includes advanced operating features and pulse timer capacity.

Easy Online Ordering at GeriatricMedical.com

404 3rd Street, Everett, MA 02149
Phone: 800-442-1205 Fax: 866-829-0506

