

SteriGear

The Fig Leaf™

Reference number 10270
Ships 20 units per carton

The Fig Leaf™ is a urinary drain bag that hides fluid from view with a built-in Fig Leaf Cover™.

The Fig Leaf™ Advantage:

- Restores patient dignity.
- Improves the environment for caregivers and visitors.
- Helps care facilities meet "quality of life" requirements.
- Saves money by eliminating extra so-called "privacy bags".
- Saves caregiver time and frustration.



FEATURES:

- 2000 ml
- Fig Leaf Cover™
- QuickView Ports™
- Drip chamber
- Non-return valve
- Anti-reflux
- Needleless sampling port
- Two air vents
- Universal hanger
- Rope hanger
- Drain valve holder
- "No Drip Flip" outlet tube
- Inlet tube
- Sampling pinch clamp
- Bed sheet clamp
- Sterile
- Latex free



Easy Online Ordering at GeriatricMedical.com

404 3rd Street, Everett, MA 02149
Phone: 800-442-1205 Fax: 866-829-0506

