



Hoyer®
Sling Guide

Information About the Caregiver(s)

When selecting slings, it is important to be mindful of the person(s) who will be fitting the sling and using the lift. Their abilities and limitations - especially experience and training - must be considered. Some slings are more complicated than others to fit and a caregiver may be reluctant to use a sling if the process of fitting it is too complex. The number of caregivers available to fit the sling as well as the total number of caregivers involved in a residents' care may also be important considerations - a range of caregivers may have a range of capabilities.

Information About Other Equipment


It is necessary to establish the sling's compatibility with other equipment, especially the lift and the lifting system. It is also important to identify any associated equipment involved in transfers, such as seating, beds, gurneys, toilets, and therapy equipment, to ascertain that the lift and sling will work effectively with these items.

Selecting the Type of Sling for the Most Common Task

To simplify the process of sling selection, **HOYER** has developed a chart to guide you in the appropriate use of each sling and task. There will, of course, be exceptions to this guide, but this reference enables a view of what tasks can be completed.

1) Use only with a net version

Sling Type	HOYER Lift	Toilet	Commode to bath ¹	Chair to Chair	Bed to Chair	Floor to Bed	Chair to Commode
Access	Advance/Presence	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Quick Fit Del	Advance/Presence	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Quick Fit	Advance/Presence	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Full Back	Advance/Presence	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Comfort	Stature	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Standing	Arise	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Transport	Arise	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>



About Joerns Healthcare

Joerns Healthcare provides the complete solution to moving and handling with a range of compatible lifts, slings and accessories designed to meet the needs of both resident and caregiver. Besides **Hoyer** slings, the Joerns Healthcare moving and handling range includes:

- Ceiling mounted lifts
- Mobile lifts
- Extended warranty and service options
- Moving and handling training packages

We trust you have found this guide both useful and informative.



Easy Online Ordering at GeriatricMedical.com

404 3rd Street, Everett, MA 02149
Phone: 800-442-1205 Fax: 866-829-0506

