MATERIAL SAFETY DATA SHEET

SECTION I — GENERAL INFORMATION

Manufacturer's Name: SMITH & NEPHEW INC - - WOUND MANAGEMENT DIVISION

Emergency Phone Number: (813) 392-1261 **Date Prepared:** DECEMBER 1996

Address: 11775 Starkey Road, Largo FL 33773-4727

Chemical Name & Synonyms: N/A

Trade Name & Synonyms: Nursing Care O Protective Ointment

Chemical Family: N/A Formula: N/A

SECTION II — HAZARDOUS INGREDIENTS / IDENTITY INFORMATION

Hazardous Components	OSHA PEL	ACGIH TLV	Other Limits	%
None	(8hr-TWA)	(8hr-TWA)		
None				

SECTION III — PHYSICAL / CHEMICAL CHARACTERISTICS

Boiling Point: > 200 Vapor Pressure (mm Hg): N/A Vapor Density (Air = 1): N/A	Specific Gravity: Melting Point: Evaporation Rate:	0.813 40–60°C < 1 (nBuAc = 1)
Solubility in Water: Negligible	pH @ 20°C:	5.6
Appearance and Odor: Plue got with a clove odor		

Appearance and Odor: Blue gel with a clove odor.

SECTION IV — FIRE & EXPLOSION HAZARD DATA

Flash Point (Method Used): $> 200^{\circ}F$ (TCC)

Flammable Limits (Lower/Upper): N/A

Extinguishing Media: Carbon dioxide, dry chemical, and foam.

Special Fighting Procedures: Firefighters should wear self-contained breathing apparatus and full protective

clothing when fighting fires.

Unusual Fire & Explosion Hazards: Use a water spray to cool exposed containers.

SECTION V — REACTIVITY DATA

	Unstable		Conditions to Avoid: Excessive heat.
Stability	Stable	Х	
Incompatibility (Materials to Avoid):):	Strong oxidizing agents.
Hazardous Decomposition products:		s:	Oxides of carbon.
Horordono	May Occur		Conditions to Avoid: See above
Hazardous Polymerization	Will Not Occur	Х	

MATERIAL SAFETY DATA SHEET

SECTION VI — HEALTH HAZARD DATA

Route(s) of Entry: Ingestion and eye contact.

Health Hazards: Irritating to eyes. Ingestion of large quantities can be hazardous.

Carcinogenicity: None

Signs and Symptoms of Exposure: Irritation of eyes and gastrointestinal tract.

Medical Conditions Generally Aggravated by Exposure: Pre-existing skin disorders may be irritated by this product.

Emergency First Aid Procedures:

EYES: Flush with large amounts of water, including under eyelids for 15 minutes. INGESTION: Give victim water and contact a physician or the poison control center for

instructions.

SECTION VII — PRECAUTIONS FOR SAFE HANDLING AND USE

Steps to be Taken in Case Material is Released or Spilled:

Wearing appropriate protective clothing, absorb spill with sand or vermiculite and place in a proper container for disposal.

Waste Disposal Method: Dispose of via a licensed waste disposal company in accordance with all local, state,

and Federal regulations.

Precautions to be Taken in

Handling and Storing:

Store in a cool, dry location away from sources of heat.

Other Precautions: N/A

SECTION VIII — CONTROL MEASURES

Respiratory Protection: N/A

(Specify Type)

Ventilatio Local Exhaust: X Special: N/A

n Mechanical: Other:

Protective Gloves: N/A Eye Protection: Safety glasses or goggles.

Other Protective Clothing or Equipment: Eyewash stations.

Work/Hygienic Practices: Normal