

# MATERIAL SAFETY DATA SHEET

PDI, THE HEALTHCARE DIVISION OF  
NICE-PAK PRODUCTS, INC.  
TWO NICE-PAK PARK  
ORANGEBURG, NY 10962  
845-365-1700

REVISION DATE: 6/03

## I. PRODUCT IDENTIFICATION

PRODUCT TRADE NAME: SANI-WIPE - ITEM# Q94384, Q73484, U01195

EPA Reg. No. 9480-7

### HAZARD RATING (NFPA)

HEALTH: 0

FLAMMABILITY: 2

REACTIVITY: 0

SPECIFIC: NONE

EMERGENCY OR INFORMATION TELEPHONE NO: 845-365-1700 (M-F DAYTIME) AT OTHER TIMES, CONTACT THE LOCAL POISON CONTROL CENTER

CHEMICAL NAME: PROPRIETARY MIXTURE

## II HAZARDOUS INGREDIENTS PER 29 CFR 1910.1200

Hazardous Ingredients	%	ACGIH TLV	CAS NUMBER
ISOPROPYL ALCOHOL	5.4800	400 ppm	67-63-0

## III-PHYSICAL/CHEMICAL CHARACTERISTICS

COLOR/ODOR/APPEARANCE: TOWELETTE SATURATED WITH A CLEAR LIQUID WITH SLIGHT ALCOHOL ODOR

BOILING POINT: 181<sup>o</sup> F (LIQUID ONLY)

VAPOR DENSITY N/A

EVAPORATION RATE: N/A

SOLUBILITY IN WATER: COMPLETE FOR LIQUID. TOWEL IS NOT SOLUBLE

SPECIFIC GRAVITY (H<sub>2</sub>O = 1): 0.990

## IV FIRE & EXPLOSION HAZARD DATA

FLASH POINT(Method Used): 112<sup>o</sup> F TAG CLOSED CUP (LIQUID ONLY) : LEL: N/A UEL: N/A

EXTINGUISHING MEDIA: WATER, FOAM

SPECIAL FIRE FIGHTING PROCEDURES: TREAT PRODUCT AS COMBUSTIBLE

UNUSUAL FIRE AND EXPLOSION HAZARDS: NONE

## V-REACTIVITY DATA

STABILITY: STABLE

INCOMPATIBILITY: ANIONIC SURFACTANT, SOAPS

HAZARDOUS DECOMPOSITION OR BYPRODUCTS: NONE KNOWN

POLYMERIZATION: WILL NOT OCCUR

CONDITIONS TO AVOID: NONE KNOWN

SANI-WIPE

**VI-HEALTH HAZARD DATA**

EFFECTS OF OVER EXPOSURE

SKIN : NONE KNOWN

EYES : Avoid Contact with eyes

INHALATION: NOT NORMAL ROUTE OF ENTRY

INGESTION: NOT NORMAL ROUTE OF ENTRY

EMERGENCY AND FIRST AID PROCEDURES:

SKIN CONTACT: If rash or irritation develops, discontinue use

EYE CONTACT: See Label Instructions

INHALATION: Remove to fresh air

INGESTION: Consult physician

TARGET ORGANS; None known

KNOWN CARCINOGEN: None known

**VII-SPILL AND DISPOSAL PROCEDURE**

SPILL CONTROL: PICK UP WITH ABSORBENT

WASTE DISPOSAL METHOD: OBEY LOCAL, STATE AND FEDERAL REGULATIONS

HANDLING AND STORAGE: AVOID HIGH HEAT AND SOURCES OF IGNITION

**VIII-CONTROL MEASURES/PROTECTION**

RESPIRATION: Not required

VENTILATION: Recommended

PROTECTIVE GLOVES: Not required with unadulterated product under normal use conditions.

EYE PROTECTION: Avoid contact with eyes.

HYGIENIC PRACTICES: Consult label for specific instructions

OTHER: Consult label for specific instructions

**IX TRANSPORT/SHIPPING**

DOT SHIPPING NAME:N/A

TECHNICAL SHIPPING NAME: N/A

DOT SHIPPING CLASSIFICATION: NOT DOT REGULATED

DOT ID NO.: N/A

DOT ;LABEL; REQUIREMENTS: N/A

UN/NA NUMBER REGULATIONS: N/A

REPORTABLE QUANTITY: N/A

**X DISCLAIMER**

THE INFORMATION FURNISHED HEREIN IS BELIEVED TO BE ACCURATE AND REPRESENTS THE BEST DATA CURRENTLY AVAILABLE TO US. NO WARRANTY, EXPRESSED OR IMPLIED IS MADE AND NICE-PAK PRODUCTS, INC. ASSUMES NO LEGAL RESPONSIBILITY OR LIABILITY RESULTING FROM ITS USE.

N/A = NOT APPLICABLE