

# ASSURE<sup>®</sup> PLATINUM

**+** PROFESSIONAL BLOOD GLUCOSE MONITORING SYSTEM



## Highly Accurate<sup>1</sup> Blood Glucose Testing in a Feature-focused System

Tight glycemic control can reduce the risk of diabetes-related complications by at least 42 percent.<sup>2</sup> The Assure<sup>®</sup> Platinum Blood Glucose Monitoring System provides healthcare professionals with highly accurate<sup>1</sup> test results, while keeping their residents' care in mind.

**arkray USA, INC.**<sup>+</sup>  
*Your Diabetes Health Ally*

# ASSURE<sup>®</sup> PLATINUM

## PROFESSIONAL BLOOD GLUCOSE MONITORING SYSTEM

The Assure Platinum Blood Glucose Monitoring System is designed for use by the healthcare professional in a multi-resident setting. A small, 0.5 µL, sample is all that is needed for each test. The LCD backlit display was developed to increase visibility in low-light situations, while the test strip ejector was engineered for touch-free test strip disposal.

### ADVANCED FEATURES

- + **qcProGuard<sup>®</sup>**: Control solution test notification provides additional quality control
- + **Hypoglycemic Alert**: Audio alert is designed to improve resident safety
- + **ProGrip<sup>®</sup>**: Rubberized-case design for secure handling

### SPECIFICATIONS

#### + **Test Strips**

Sample size: 0.5 µL  
Test time: 7 seconds  
Result range: 20-600 mg/dL  
Chemistry: Glucose oxidase

#### + **Technology**

Calibration: Plasma-referenced  
Auto Coding: Yes  
Display: LCD  
Memory: 500 tests with time and date stamp

#### + **Meter**

Size: 4.5 x 2.5 x 1.2 Inch  
Weight: 2.8 oz. without batteries

#### + **Power**

Battery: Two AAA batteries  
Tests/battery: More than 5,000



### COMMITMENT TO QUALITY

ARKRAY offers high quality products and programs to better manage the disease for people living with diabetes and their healthcare providers. The Assure portfolio is designed to deliver optimal quality while remaining affordable in today's environment.

### ORDER DETAILS

Product Name	Quantity	Product Number	UPC Number 015482-	NDC Number 08317-
Assure Platinum Blood Glucose Meter	1	500001	50001-6	5000-01
Assure Platinum Test Strips	50	500050	50050-4	5000-50
Assure Platinum Test Strips	100	500100	50100-6	5001-00
Assure Dose Control Solution – Level 1 (N) and 2 (H)	1	500006	50006-1	5000-06

#### References

1. ARKRAY, Inc. Assure Platinum 510(k) substantial equivalence determination decision summary assay and instrument combination template. Dec 2009;1-10.
2. DCCT and EDIC: The Diabetes Control and Complications Trial and Follow-up Study <https://www.niddk.nih.gov/about-niddk/research-areas/diabetes/dcct-edic-diabetes-control-complications-trial-follow-up-study/Pages/default.aspx#eye>, NIH Publication No 08-3874, May 2008

ARKRAY USA, INC. <sup>+</sup>

Minneapolis, MN 55439 • 800.818.8877 • [www.arkrayusa.com](http://www.arkrayusa.com)

©ARKRAY USA, Inc. All rights reserved.  
AT006-01 Rev. 11/17