

# ASSURE<sup>®</sup> PRISM MULTI

+ BLOOD GLUCOSE MONITORING SYSTEM



## Fast and Accurate<sup>1</sup> Blood Glucose Testing in an Easy-to-use Meter

One in every three residents in long-term care have diabetes.<sup>2</sup> Nationally, it costs an additional 14 billion dollars each year to care for these residents compared to residents without diabetes.<sup>3</sup> The Assure<sup>®</sup> Prism multi Blood Glucose Monitoring System provides healthcare professionals with a fast, accurate and easy way to monitor their residents' blood glucose levels.

ARKRAY USA, INC.<sup>+</sup>  
*Your Diabetes Health Ally*



# ASSURE<sup>®</sup> PRISM MULTI

## BLOOD GLUCOSE MONITORING SYSTEM



### RELIABLE DIABETES CARE

The Assure Prism multi Blood Glucose Monitoring System has a clinically driven design and advanced technology. This multi-resident meter complies with ISO 15197:2013 accuracy requirements.<sup>4</sup> Design features like a test strip that points toward the resident and backlit display are designed to make it easier for the healthcare professional to monitor a resident's blood glucose level.

### PERFORMANCE FOCUS

- + Test Results in 5 Seconds**  
Fast test results reduce downtime
- + Optimized Readout Display**  
Large, thick numbers are intended to improve readability
- + Compact Design**  
Small, lightweight design provides better portability
- + Test Strip Release Button**  
Safely dispose of used test strips

### SPECIFICATIONS

#### Test Strips

Sample size: 0.5 µL  
 Test time: 5 seconds  
 Result range: 20-600 mg/dL  
 Chemistry: Glucose oxidase

#### Technology

Calibration: Plasma-equivalent  
 Auto Coding: Yes  
 Display: LCD  
 Memory: 500 tests

#### Meter

Size: 3.66 x 2.20 x 0.83 Inch  
 Weight: 2.1 oz. with batteries

#### Power

Battery: Two 3.0V Lithium batteries  
 Tests/battery: 1,000

### ORDER DETAILS

Product Name	Quantity	Product Number	UPC Number 015482	NDC Number 08317
Assure Prism multi Blood Glucose Meter	1	530001	53001-3	5300-01
Assure Prism multi Test Strips	50	530050	53050-1	5300-50
Assure Prism multi Test Strips	100	530100	53100-3	5301-01
Assure Prism Control Solution – Level 1 and 2	1	530006	53006-8	5300-06

#### References

1. ARKRAY Assure Prism multi Blood Glucose Monitoring System. User Instruction Manual. 5118-03 Jan 2015. 45-7.
2. Boschert S. Diabetes found in one-third of residents. Caring for the Ages. 2008;9(1):23.
3. Yang W, Dall T, Halder P, et al. Economic cost of diabetes in the U.S. in 2012. ADA. 2013;36(1):1033-46.
4. Walker J, Gill P, Gleisner J. Accuracy of the Assure Prism multi as it relates to the ISO 15197:2013 requirements in the monitoring of diabetes mellitus. 11th World Congress on Insulin Resistance, Diabetes and Cardiovascular Disease. 2015. Los Angeles, United States. AS003-00.

ARKRAY USA, INC.

Minneapolis, MN 55439 • 800.818.8877 • www.arkrayusa.com

© ARKRAY USA, Inc. All rights reserved.  
 AS006-02 Rev. 09/17