

Biatain[®] Silicone Lite
Foam Dressing

Enjoy the freedom of the
new Biatain[®] Silicone Lite



Biatain® Silicone Lite

– superior absorption with soft adhesion for increased mobility

Features

- Unique Biatain foam structure for superior absorption
- Secure fit
- Superior exudate management¹
- Ultra thin, soft and flexible
- Aseptic 3-piece non-touch opening




Benefits

- Conforms closely to the wound bed and body
- Increased mobility
- Minimal maceration and leakage
- Easy-to-use application reduces wasted dressings

The new and improved Biatain Silicone Lite

- The new silicone adhesive wound contact layer provides a secure fit to keep the dressing in place while still ensuring minimal pain upon removal.
- It is a thin, flexible foam dressing that can be used when there is a need for increased mobility.
- Biatain Silicone Lite has a 3-piece non-touch application.




Usage

Exudate Level

Biatain Silicone Lite is indicated for a wide range of acute and chronic exuding wounds such as leg ulcers, pressure ulcers, non-infected diabetic foot ulcers, donor sites, post-operative wounds and skin abrasions. Biatain Silicone Lite is suitable on low-to-moderately exuding wounds and can also be used in combination with compression therapy.

New Biatain Silicone products



Biatain Silicone

– superior absorption with soft adhesion for general purposes



Biatain Silicone Lite

– superior absorption with soft adhesion for increased mobility

Biatain Silicone Lite

Item code	Size in inches	Units per box	HCPCS*
33452	2x2"	5	A6212
33444	3x3"	10	A6212
33445	4x4"	10	A6212
33446	5x5"	10	A6212
33453	2x5"	5	A6212

1. Andersen MB, Marburger M. Comparison of absorption and fluid handling between ten wound dressings with silicone adhesive. Wounds UK. Harrogate, UK 2014.

*Reimbursement Disclaimer: Coloplast Corp. provides this information for your general reference and related to the reimbursement of Coloplast products only. Reimbursement, coverage and payment policies can vary from one insurer and region to another, and may change over time. Coloplast does not guarantee coverage or payment of products.