

# GLUCOCARD<sup>®</sup> EXPRESSION<sup>™</sup>

+ BLOOD GLUCOSE MONITORING SYSTEM



See the Quality.  
Hear the Results.

With our advanced, yet aggressively priced blood glucose meters, you'll worry less about costs so you can cover more customers while maintaining your profits.

Leverage our products and programs to help your customers better control their diabetes and live richer, healthier lives. Become the #1 destination for their diabetes care.

arkray USA, INC. +  
*Your Diabetes Health Ally*



# GLUCOCARD<sup>®</sup> EXPRESSION<sup>™</sup>

BLOOD GLUCOSE MONITORING SYSTEM

GLUCOCARD Expression is an audio-enabled meter with a large display and tactile buttons ideal for those with dexterity or visual impairment. GLUCOCARD Expression can help people with diabetes see how making small, simple lifestyle changes can have a big impact on their health.

## SPECIFICATIONS

### + *Test Strips*

Name: GLUCOCARD Expression

Sample size: 0.8 µL

Test time: 6 seconds

Meter range: 20-600 mg/dL

Chemistry: Glucose oxidase

### + *Technology*

Results: Plasma referenced

Calibration: Auto code

Audio: English and Spanish

Memory: 300 tests with time and date stamp

Averaging: 7-, 14- and 30-day averaging

### + *Meter*

Size: 3.7 x 2.3 x 0.8 inches

Weight: 2.8 oz

Battery: 2 AAA batteries

Tests/Battery: 1,000 tests of continuous use or 1 year



## ORDER DETAILS

Product Name	Count	UPC Number	NDC Number
GLUCOCARD Expression Meter Kit	1	015482-57110-8	08317-5711-00
GLUCOCARD Expression Meter Basic	1	015482-57010-1	08317-5700-01
GLUCOCARD Expression Test Strips Retail	50	015482-57205-1	08317-5720-50
GLUCOCARD Expression Test Strips DME	50	015482-57005-7	08317-5700-50

ARKRAY USA, INC.<sup>+</sup>

Minneapolis, MN 55439 • 800.566.8558 • www.arkrayusa.com

©ARKRAY USA, Inc. All rights reserved.  
GX002-02 Rev. 08/17