

CHECK OUT OUR NEW & IMPROVED

GRAVITY SERIES PRESSURE REDISTRIBUTION MATTRESSES

THERAPEUTIC FOAM MATTRESSES WITH IMPROVED
QUALITY AND FEATURES FOR THE SAME PRICE!


Gravity 7 Long Term Care Pressure Redistribution Mattress

- Therapeutic die cut foam now extends across the entire 36" width creating a greater patient contact surface area. Firm high density perimeters remain beneath for safe patient transfers.
- Improved foam densities increase weight capacity without increasing interface pressure. Weight cap now 350 lbs!
- CFR 1632/CFR 1633 Fire rating coded
- Available in 76", 80" and 84" Lengths
- Available with or without raised side rails
- Warranty: 7 years foam / 1 year cover

Item

15876:	36" (W) x 76" (L) x 6" (H)	15877:	36" (W) x 76" (L) x 6" (H) W/ RAISED RAILS
15770:	36" (W) x 80" (L) x 6" (H)	15777:	36" (W) x 80" (L) x 6" (H) W/ RAISED RAILS
15784:	36" (W) x 84" (L) x 6" (H)	15785:	36" (W) x 84" (L) x 6" (H) W/ RAISED RAILS

350 LB
Weight Capacity




Gravity 8 Deluxe Long Term Care Pressure Redistribution Mattress

- Therapeutic die cut foam now extends across the entire 36" width creating a greater patient contact surface area. Firm high density perimeters remain beneath for safe patient transfers.
- Improved foam densities increase weight capacity without increasing interface pressure. Weight cap now 400 lbs!
- CFR 1632/CFR 1633 Fire rated
- Available in 76", 80" and 84" Lengths
- Available with or without raised side rails
- Warranty: 10 years foam / 1 year cover

Item

15886:	36" (W) X 76" (L) X 6" (H)	15888:	36" (W) X 80" (L) X 6" (H) W/ RAISED RAILS
15870:	36" (W) X 80" (L) X 6" (H)	15885:	36" (W) X 84" (L) X 6" (H) W/ RAISED RAILS
15884:	36" (W) X 84" (L) X 6" (H)		

400 LB
Weight Capacity




Gravity 9 Premium Long Term Care Pressure Redistribution Mattress

- Therapeutic visco elastic memory foam now extends across the entire 36" width creating a greater patient contact surface area. Firm high density perimeters remain beneath for safe patient transfers.
- Improved foam densities increase weight capacity without increasing interface pressure. Weight cap now 450 lbs!
- CFR 1632/CRF 1633 Fire rated
- Available in 76", 80" and 84" Lengths. COMING SOON 42" WIDTH!
- Available with or without raised side rails
- Warranty: 15 years foam / 1 year cover

Item

15996:	36" (W) X 76" (L) X 6" (H)	15997:	36" (W) X 76" (L) X 6" (H) W/ RAISED RAILS
15970:	36" (W) X 80" (L) X 6" (H)	15977:	36" (W) X 80" (L) X 6" (H) W/ RAISED RAILS
15984:	36" (W) X 84" (L) X 6" (H)	15985:	36" (W) X 84" (L) X 6" (H) W/ RAISED RAILS

450 LB
Weight Capacity



CUSTOM SIZES AVAILABLE

