



# HD2 Bladder Scanner

*World's Most Affordable  
Bladder Scanner*



With the **HD2** BLADDER SCANNER, we bring you:

### ACCURACY

Patented Algorithm to guarantee  
 $\pm 12\%$  clinical precision



*Reliable in the measurement of bladder volume and post void residual. Idea solution for nursing and long-term care facilities.*

### AUTOMATIC VOLUME CALCULATION

One Touch, Calculate All



*Efficient Solution to your regular nursing work. Search "HD2" in Youtube to watch the easy training video!*

### NO NEED FOR CALIBRATION

in the Entire lifespan by  
HD2 Technology



*We take your cost saving into consideration not only in the beginning, but for all the time.*

### TWO PRE-SCAN MODES

Easy Mode  
Expert Mode



*We respect your habit of bladder scanning in the mean time providing abundant clinical information.*

### VOICE RECORDING FUNCTION URINE VOLUME SETTING AND ALARM SETTING



*Spare your hand eliminate sophisticated procedures; Intuitively showing the catheterization timeline.*

### PORTABLE

3.5 hours  
for continuous scanning

Mode	Expert
Sound	Large >
Remind	OFF >
Language	English >
Calibration	OFF >
Update	>
Version	HD2_V1.1.0.0

*8 hours for stand-by; SD card storage; Longer time using; Unlimited saving capacity for patient cases.*



CALL FOR MORE INFO  
**888.602.PROS**



The HD2 Scanner aims to improve healthcare and create value in the medical field by researching, manufacturing and marketing our bladder scanners product line globally



The medical cost of Healthcare-associated infections (HAIs) in the United States is as high as \$45 billion per year, causing big social problems. (According to the data of Centers for Disease Control and Prevention, CDC in 2014)

As the most common type of medical infection, catheter-associated Urinary Tract Infections account for up to 40% of HAIs. This high incidence of CAUTIs causes pain and severe clinical risk in the patients.

Super-Affordable, Portable **HD2 Bladder Scanner** provides quick and accurate measurement for bladder volume and post-void residual.

With the **HD2 Bladder Scanner**, we bring the value

- ▶ Helps prevent unnecessary catheterization, reducing catheter associated urinary tract infections
- ▶ Helps reduce the incidence of UTI-related hospital stays
- ▶ Saves the UTI-associated medical cost
- ▶ Improves the welfare and satisfaction of the patients and aged-care customers.

## With HD2, We Create Comfort and Value Less Pain for Patients | Easier Work For Nurses

SCANNING MODE	Mechanical Sector-shaped Scanning
FREQUENCY OF PROBE	3.5 MHz
VOLUME MEASUREMENT RANGE	0-999ml
SIZE OF MONITOR	2.5 inch TFT LCD
MODE	Easy Mode
OTHER FUNCTION	Voice Recording, Alarm function
LANGUAGE	Multi-Language
BATTERY CAPACITY	2600mAh
STORAGE	TF Card information storage
DIMENSIONS	203*213*55mm
WEIGHT	approximate 500g
CONTINUOUS SCANNING TIME OF BATTERY	>3.5 hours
STAND BY TIME OF BATTERY	>8 hours

  
HD2 BLADDER SCANNER

CALL FOR MORE INFO  
**888.602.PROS**