



**To Whom it May Concern:**

**Stander's Bed Rail Advantage, item #5000, protects against the seven entrapment zones outlined by the FDA, when properly installed. Proper installation includes:**

- **Safety Strap secured tightly to bed frame so that the mattress is unable to slide away from the Bed Rail Advantage.**
- **Placement of bed rail should be no closer than 18" to headboard or footboard**
- **Pouch Organizer is used over the handle to prevent entrapment. The organizer snaps securely into place.**
- **The height setting of the Bed Rail Advantage must be adjusted such that there is no gap between the top of the mattress and the bottom of the handle.**
- **Two Bed Rail Advantages must not be placed side-by-side on the same side of the bed.**

**Please feel free to contact us if you have any further questions of if we can assist you in any way.**

**Sincerely,**

**Brandon Birch  
COO  
Stander, Inc**